



PRECISION BALANCE **BJA**

User Manual



v.201811




HiWEIGH
Weighing system & solution

An ISO9001 registered company
@No.335 Haishen, Xingxin Road, Huinan Town, Pudong District, Shanghai 201301, China
www.hiweigh.com All rights reserved, specifications subject to change without notice

Value Each Gram

Before Use

Surroundings required for operating the balance




1. Avoid corrosion, the room temperature should be stable.
2. Place on a stable desk, professional balance table is recommended.
3. Place in a dry atmosphere, avoid strong air flow and intensive magnetic field.

Preparing for use

1. Unpack the package, take out the balance.
2. Clean the appliance with soft cloth or brush.
3. Put the s.s pan on the balance.
4. Connect the balance to the power supply.


Preheating

In order to obtain the accuracy of its use, the balance should be preheated for no less than 30 minutes after it is connected to the power.




Keypad and Display

Parts and Buttons



1. DISPLAY 2. ON/OFF 3. OUTPUT 4. PCS 5. PER (PERCENTAGE) 6. CAL 7. TARE 8. PAN

Display



1. OPERATING INDICATOR 2. +,- 3. WEIGHT 4. UNIT

MAIN PARAMETER AND SURROUNDINGS

Technical Parameter

Model	BJA103	BJA123	BJA203	BJA303	BJA403	BJA503
Capacity(g)	110	120	210	310	410	510
Real Scale Division(g)	0.001	0.001	0.001	0.001	0.001	0.001
Examined Scale Division	0.01	0.01	0.01	0.01	0.01	0.01
Standard Weight(g)	100	100	200	300	400	500
Output Interface	RS232					
Main Function	Tare, calibration, percentage, unit exchange, counting, automatic malfunction check.					
Dimensions	320×200×320mm					
Power	AC110V50Hz or AC220V50Hz 14VA					
Surroundings	10°C-30°C ≤85%RH					

Model	BJA502	BJA1102	BJA2102	BJA3102
Capacity(g)	510	1100	2100	3100
Real Scale Division(g)	0.01	0.01	0.01	0.01
Examined Scale Division	0.1	0.1	0.1	0.1
Standard Weight(g)	500	1000	2000	3000
Output Interface	RS232			
Main Function	Tare, calibration, percentage, unit exchange, counting, automatic malfunction check.			
Dimensions	320×200×320mm			
Power	AC110V50Hz or AC220V50Hz 14VA			
Surroundings	10°C-30°C ≤85%RH			

OPERATIONS

Calibration

After preheating, before use, it is necessary to calibrate the balance to obtain the weighting accuracy.

1. Press TARE

Remove the items on the pan, press TARE to make sure the display of the balance is zero (E.g.: BJA502).

2. Press CAL

Press CAL, the display will be "C".

3. Put the standard weight

Put the standard weight on the pan according to Technical Parameter.

4. Finish Calibration

When the display of the balance equal to the weight of the standard weight and it comes with a sound of du, which indicates the calibration is success, the balance will return to the weighing mode automatically. Then you can remove the weight. (E.g.: BJA502)

Caution: After pressing CAL, if "CE", an error sign, displays, it indicates the balance is unstable or there are some items on the pan, then remove the items and then press TARE for re-calibrating.

Weighing

1. Put required weighting goods on the pan and wait till the stable sign "O" appears, then you can read the weight of the goods.
2. If second item is required, press TARE to zero the balance.
3. Put the second item on the pan, the display will be the weight of the 2nd item.
4. Press TARE again to make sure the display is zero.
5. Remove the two items on the pan, the display will be total weight of the two items.

Counting

1. Remove the items on the pan, press TARE to make sure the display is zero.
2. Put 10 pieces of the required items on the pan, press PCS after the reading is stable, the balance will turn into the counting mode, "10" will be displayed on the balance.
3. Add or remove several items on the pan, the quantity will display automatically.
4. If you want to know the weight of the items on the pan, just press UNIT, it will exchange the quantity to the weight value. Press UNIT again, it will back to quantity display again.

5. If you want to exit the counting mode, just press PCS again in that mode.

Cautions: 1). *Different types of items should be set respectively.*
2). *Each weight of counting items should be more than 10d, and the difference of the weight of the items should be less than the accuracy of the balance.*
3). *The ex-factory sample size for the counting function is set to 10, and there are also 20,40,80, three shifts for selecting.*

Percentage

The weight ratio between the two weights will be calculated on the basis of the percentage counting function of the balance.

1. Remove the items on the pan, press TARE to make sure the display is zero
2. Put the reference sample on the pan, after the reading is stable, press PERCENTAGE, the balance is in the PERCENTAGE mode, then the display will be "100%".
3. Add or remove the goods on the pan, the data displayed on the balance is the percentage of the weight of the goods on the pan and that of the reference sample
4. If you want to know the percentage of the weight of the added or removed goods and that of the reference sample, press TARE first to make sure the display is "0%", then add or remove the goods.
5. If you want to know the weight of the goods on the pan, press UNIT, then the percentage will be changed into weight display, then press UNIT again, the weight will be changed into percentage.
6. If you want to exit the PERCENTAGE mode, just press PERCENTAGE again in that mode.

Cautions 1). *Different types of items should be set respectively.*
2). *The weight of the reference sample should be no less than 30% of the permitted maximum weight, otherwise, "E1" will be displayed, in that case the sample should be reset.*

Unit

The conversion function of this balance is available, such as "g" can be changed into "ct".

1. Remove the items on the pan, press TARE to make sure the display is zero.
2. Put one required item on the pan, when "g" displays, it indicates the reading is stable, then the display of the balance is the weight of this item.
3. Press UNIT, the display will be the value of "ct"
4. Then press UNIT again, unit "ct" will be changed into "g".

Output interface

With RS-232 output interface, the balance can be connected to other instruments

such as computer, electronic balance recording instrument.

Trouble Shooting

Symptoms	Possible cause	Solution
No display	power disconnected.	Check the power cord.
The display of L	The pan is placed improperly	Re-place the pan
	No calibrating for the balance.	Re-calibrate the balance.
The display of H	Overload	Remove the items on the pan.
	Error for balance calibration.	Re-calibrate the balance.
No stable display	The pan is placed improperly	Re-place the pan
	Vibration or airflow.	Choose proper surroundings.
Big error of display	No tare before weighing.	Tare before use.
	No calibrating for the balance.	Re-calibrate the balance.

Full set of products

- | | | |
|-------------------------|---|---------------------------------|
| 1. Main unit of balance | 1 | |
| 2. Pan | 1 | |
| 3. Power cord | 1 | |
| 4. Instruction Manual | 1 | |
| 5. Standard weight 100g | 1 | <i>(Only for BJA103 BJA123)</i> |
| 6. Standard weight 200g | 1 | <i>(Only for BJA203)</i> |



v.201811
User Manual
BJA